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What positioning and fall alarms can and cannot promise

Positioning and automatic fall detection can provide important information, but they are not absolute guarantees. Coverage, environment, wearing position, device status and configuration affect what the recipient can actually see.

PUBLISHED	12 September 2026
AUTHOR	Siren Tech's engineering team
READING TIME	6 min read

Risk scenario: the alarm feature exists but the chain is untested

Imagine a user wearing a watch with SOS and fall detection. The recipient has app access, but no one has tested indoor positioning, who responds outside normal hours or what happens if an automatic fall alarm is not triggered.

This is a risk scenario. It shows why a feature list is not proof that the complete safety workflow works in its intended environment.

Test where the device will be used

Position can be affected by coverage and whether the user is indoors or outdoors. Fall detection is affected by how the device is worn and the type of movement. Test SOS, calls, position and automatic alarms in the real environments before normal operation.

Also decide how low battery, lost contact or an uncertain position should be handled by the recipient.

Four decisions before launch

Document responsibilities and limitations together with users and recipients.

- Who receives SOS, fall and inactivity alarms?
- Which position and history may each role see?
- What fallback procedure applies when coverage or battery is unavailable?
- How often are the device and recipient chain tested?

How Siren Appen can contribute

Siren Appen can bring configured SOS, fall and inactivity alarms, available position, geofences, history, calls and device status together for authorised recipients. Exact functionality depends on the compatible device and agreed workflow.

Supported health values are indicative and do not replace medical equipment, assessment or diagnosis. Personal data and location access require a clear purpose and appropriate permissions.

Sources and further reading

[Book a Siren Appen demo](#)

[See compatible personal alarms](#)

[Swedish Authority for Privacy Protection: data protection principles](#)